



Silk'n[®]
Beautiful Technology

FaceTite[™]

USER MANUAL

Home Skinovations Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Home Skinovations Ltd. is believed to be accurate and reliable at the time of publication. However, Home Skinovations Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Home Skinovations Ltd. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Home Skinovations Ltd. Data is subject to change without notification. Home Skinovations Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Home Skinovations Ltd. Specifications are subject to change without notice.

Home Skinovations Ltd. Tavor Building, Shaar Yokneam POB 533, Yokneam 2069206, ISRAEL

www.silkkn.co.il, contact@silkkn.co.il

Table of Contents

1. Before you start	2
1.1. Silk'n FaceTite device description	2
1.2. Package contents	2
2. Silk'n FaceTite Intended Use	3
3. Safety with Silk'n FaceTite.....	3
4. Contraindications	4
5. Do not use Silk'n FaceTite if:	5
6. What to expect with Silk'n FaceTite.....	5
7. Possible side effects.....	5
8. Get started!.....	6
9. Treatment regime	7
9.1 Treatment zones.....	8
10. Maintenance of Silk'n FaceTite	9
10.1 Cleaning Silk'n FaceTite.....	9
11. Troubleshooting	9
12. Customer service	10
13. Specifications	10
14. Labeling.....	11

1. Before You Start

Before using Silk'n FaceTite for the first time, please read this Instructions for Use booklet in its entirety, and watch the Silk'n FaceTite 'How to Use' video.

Pay particular attention to sections on device use procedures, device operation, and after-use procedures. We recommend you re-familiarize yourself with this Instructions for Use booklet before each use of Silk'n FaceTite.

Silk'n FaceTite is a powerful electrical device. As such, it should be used with special attention to safety. Please read all warnings and safety precautions before use, and strictly follow them when using Silk'n FaceTite.




Keep away from children!

1.1 Silk'n FaceTite device description

Silk'n FaceTite is an anti-aging face treatment device which utilizes RT™ (Home Tightening) technology. This technology combines three kinds of energy: Infrared Radiation (IR) that affects the upper layer of skin for improved skin texture, Red Light (LED) that reaches into the dermis for rejuvenation and collagen restoration, and focused Bi-Polar Radio Frequency energy (Bi-Polar RF) that penetrates deep into the subdermal fat for deep tissue tightening and wrinkle reduction.

Silk'n® FaceTite has a TREATMENT SURFACE (1) with three ELECTRODES (2) and a CENTRAL BUTTON (3) with five surrounding INDICATOR LIGHTS (4).



- ① TREATMENT SURFACE
- ② ELECTRODES
- ③ CENTRAL BUTTON
- ④ INDICATOR LIGHTS

1.2. Package contents

Upon opening the Silk'n FaceTite package, you will find the following parts:

- Silk'n FaceTite Device
- A tube of Slider Gel
- A power adaptor
- This user manual
- A beauty case

2. Silk'n FaceTite Intended Use

Silk'n FaceTite is intended for skin texture improvement and skin tightening, collagen restoration and wrinkle reduction. Silk'n FaceTite is suitable for use on your cheeks, forehead, chin, neck and the area around your eyes.

Please note: Silk'n FaceTite is intended both for self-treatment, and for the treatment of others.

3. Safety with Silk'n FaceTite

Silk'n FaceTite protects your skin

Silk'n FaceTite comes with a unique built-in heat stabilizer designed to measure skin temperature and stabilize the emitted heat energy.

To prevent overheating, a unique built-in heat stabilizer inside the Silk'n FaceTite measures the treated skin temperature during session. When a rise in the treated skin temperature is detected (42°C-43°C), the stabilizer will automatically stop emitting heat, until skin temperature lowers enough for the stabilizer to resume heating.

Note: Make sure the device is tightly pressed to your skin before starting and during treatment!

4. Contraindications

Important Safety Information – Read Before Use!

Silk'n FaceTite is not designed for everyone. Please read and consider the information in the following section before use. For further information and personalized advice you may also visit our website www.silkn.co.il.

- Do not use Silk'n FaceTite if you have a pacemaker or internal defibrillator, or any other active electrical implant anywhere in your body.
- Do not use Silk'n FaceTite if you have a permanent implant in the treated area such as metal plates and screws, silicone implants or an injected chemical substance.
- Do not use Silk'n FaceTite if you have a current history or condition of cancer, or pre-malignant moles.
- Do not use Silk'n FaceTite if you have received radiation therapy or chemotherapy treatments within the past 3 months.
- Do not use Silk'n FaceTite if you are pregnant or nursing.
- Do not use Silk'n FaceTite if you have a history of diseases stimulated by heat, such as recurrent Herpes Simplex in the treatment area. If so you may treat only following a prophylactic regime.
- Do not use Silk'n FaceTite if you suffer from poorly controlled endocrine disorders, such as diabetes and thyroid conditions.
- Do not use Silk'n FaceTite if you suffer from epilepsy.
- Do not use Silk'n FaceTite if you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria, lupus, etc.
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.
- If your skin was treated with a physician-prescribed drug within the past 6 months, please consult with your physician.
- Do not use Silk'n FaceTite if you have any active condition in the treatment area, such as sores, psoriasis, eczema, and rash.
- Do not use Silk'n FaceTite if you had any surgery in the treated area within 3 months prior to treatment.

5. Do not use Silk'n FaceTite if:

- the attachments or accessories are not recommended by Home Skinovations Ltd.
- it has a damaged cord or plug.
- you see or smell smoke when it is in use.
- it is not working properly or if it appears damaged.
- the THERMAL SENSOR is cracked, or broken.
- the outer shell is cracked or is coming apart.
- you are bathing, or if it becomes damp or wet. Unplug immediately if it has fallen into water.

Warning: Keep Silk'n FaceTite away from water! Do not place or store Silk'n FaceTite where it can fall, be pushed or placed into a tub, sink or any other vessel containing water.

6. What to expect with Silk'n FaceTite

For many people, using Silk'n Titan may be their first experience with a heat-based device designed for home use. During a treatment session it is normal to experience or feel:

- a sensation of warmth: This is from the bi-polar RF energy.
- some redness or slight swelling: This should subside within a few hours, or up to one day.

7. Possible side effects

Excessive redness or swelling: This is rare and should normally disappear in within a few days. Cool the affected area with cold water for 10 minutes, and apply a remedy against burns.

Consult a skin physician before use if you have couperose or rosecrea.

Consult your physician if you experience any persistent symptoms including pain, skin redness, swelling, or a burning sensation.

8. Get Started!

1. Remove Silk'n FaceTite DEVICE and other components from box.
2. Plug the POWER ADAPTOR cord into the Silk'n FaceTite socket.
3. Plug the POWER ADAPTOR into an electrical outlet. All five indicator lights will cycle repeatedly one after the other.
4. The device is now in STANDBY mode and is ready to start.
5. Press the CENTRAL BUTTON. The device will turn on to energy level 1 automatically. One INDICATOR LIGHT will turn on.

Setting Energy Levels:

- Silk'n FaceTite offers 5 levels of energy represented by 5 ENERGY LEVEL INDICATOR LIGHTS.
 - The selected energy level determines the intensity of Silk'n FaceTite energy delivered to your skin.
 - To set an energy level, press the CENTRAL BUTTON repeatedly to move between energies, from 1 to 5 and back to 1.
 - When the desired energy is set, the ENERGY LEVEL INDICATOR LIGHTS will indicate the level chosen.
 - To put FaceTite™ back in standby mode, hold the CENTRAL BUTTON until all five lights will cycle continuously.
 - If you wish to change energy levels mid-treatment, press the button again. This will change the energy level on the spot and allow you to treat the same or different areas with different energy levels according to your preferences.
6. After setting the energy level for the treatment, spread a layer of Silk'n® Slider Gel over the area you intend to treat.

Note: Slider Gel (or any other water-based clear gel) facilitates smooth movement of the treatment surface over the skin and enables conductivity. To ensure safe and efficient treatment, always start by applying Slider Gel!

7. Press the device to the designated skin area. You will feel a mild sensation of warmth. Move the device in gentle massage- like linear motions on the area you wish to treat. Progress slowly and carefully.

Please Note: In order for the device to work properly, it must be firmly connected to the skin!

8. Immediately after treatment apply a facial lotion, moisturizer or mask.

9. Treatment regime


First Treatment

Find the right energy level for each facial site:

- Start at Energy level 1 for a few minutes (see treatment instructions below).
- Raise the energy level gradually for a few minutes at each level.
- Stop at the highest energy level that you are still comfortable with.
- It is advisable to reach the highest energy for maximum clinical effect.

Further Treatments:

- Turn on and set to the last comfortable energy level setting.
- Apply a thin layer of gel to the treatment facial area.
- Repeat twice a week for 10 Weeks. For maintenance, treat once every 2 weeks.
- For maximum results - each treatment zone should be treated for at least 15 minutes per treatment.



9.1 Treatment Zones

1. FOREHEAD
including in-between brows
Tightening

2. RIGHT CHEEK
including around the eyes
Treat fine lines and Volume building

3. LEFT CHEEK
including around the eyes
Treat fine lines and Volume building

4. UPPER LIP
Including nasolabial fold
Tightening and Eliminating wrinkles

5. CHIN
Including jawline
Tightening and V-shaping

6. NECK
Smoothing

FOREHEAD

Tightening

CHEEKS

Volume Building

UPPER LIP

Eliminating Wrinkles

JAWLINE

V Shapping

CHIN

Tightening

AROUND THE EYES


Treat Fine Lines

NASOLABIAL FOLD

Tightening

NECK

Smoothing



10. Maintenance of Silk'n FaceTite

Note: Always keep Silk'n FaceTite stored in the beauty case that was delivered with it. Store Silk'n FaceTite and the POWER ADAPTER in a cool, dry place, out of reach from children

10.1. Cleaning Silk'n FaceTite

- Unplug Silk'n FaceTite before cleaning. After every treatment session, it is necessary to clean your Silk'n FaceTite device, and especially the TREATMENT SURFACE.
- Use absorbent paper (such as toilet paper or hand drying napkins) to gently wipe the gel residue off the TREATMENT SURFACE.
- Never immerse Silk'n FaceTite or any of its parts in water!

11. Troubleshooting

“My Silk'n FaceTite does not start.”

- Make sure the POWER ADAPTOR is properly connected to the Silk'n FaceTite device.
- Make sure the POWER ADAPTOR is plugged into an electrical outlet on the wall.

“Heat energy is not emitted and I don't feel any heat emitted from the device even though I pressed the CENTRAL BUTTON”

- Heat will not be emitted if the TREATMENT SURFACE is not connected to your skin. This is a safety feature. Make sure the device is tightly pressed to your skin before starting treatment.
- Heat will not be emitted from the TREATMENT SURFACE when not using Slider Gel, or previously applied Slider Gel has dried up. Apply Slider Gel to the area you intend to treat directly before treatment.
- Do not attempt to open or repair your Silk'n FaceTite device. Only authorized Silk'n® personnel are permitted and able to perform repairs.
- Trying to open Silk'n FaceTite will void your warranty and may also damage the device. Please contact Silk'n® Customer Service if you have a broken or damaged device in need of repair.



When Silk'n FaceTite reaches its End of Life stage, the device is to be orderly disposed of. Deliver it to the nearest electronic waste collection point for safe and efficient disposal.

12. Customer Service

More information about Silk'n products, please visit our website www.silkn.co.il.

If your device is broken, damaged, in need of repair or any other assistance, please contact our Customer Service center via our:





- Website: www.silkn.co.il
- Email: contact@silkn.co.il
- Phone: 1-700-700-309

Our customer service representatives will be happy to help you!

13. Specifications

RT™		Technology
620nm±20nm;		Wavelength Red light
850nm±20nm;		Wavelength IR
1 MHz ±30%; 12W max	Radio Frequency	
Heat stabilizer preventing skin overheating		Operation and Safety
22cm X 22cm X 9cm		Package Size
125gr		System weight
10°C to 35°C	Operating	Temperature
-40 to +70°C	Storage	
30 to 75%rH	Operating	Relative Humidity
10 to 90%rH	Storage	
700 to 1060hPa	Operating	Atmospheric pressure
500 to 1060hPa	Storage	
KSA24A1200150HE	Adaptor Model	
100-240V; 50-60Hz; 0.5A	Rated	
12Vdc; 1.5A	Output	

14. Labeling

-  Class II equipment
-  WEEE - Waste Electrical and Electronic Equipment
-  CE Mark
-  Do not use this appliance in a wet surrounding

Degree of protection against ingress water: ordinary

This device is not suitable for use in the presence of flammable anesthetic mixture with air or with Oxygen or Nitrous Oxide.



The logo features a dark grey, wavy, brushstroke-like graphic element above the text.

HomeSkinovations™

contact@silkn.co.il

www.silkn.co.il